

Experiment # 2: Part-A

Measure velocity profile, both mean and fluctuating component, in the turbulent wake behind a circular cylinder using hot-wire anemometer.

Learning Objectives:

- a) Use of hot-wire anemometry for velocity measurement
- b) Analysis of hot-wire data to obtain mean and fluctuating part of the velocities
- c) Wake characteristics behind a circular cylinder

Proposed Plan:

- a) Study the experimental set-up including each and every component with specification.
- b) Study the circuit diagram for a constant temperature hot wire circuit.
- c) Mount the hot-wire probe in the traversing mechanism and make the necessary signal connection.
- d) Make a block diagram of the experimental set up.
- e) Calibrate hot wire probe in the free stream without any model in the tunnel.
- f) Use a rough cylinder model to generate turbulent wake
- g) At one free stream wind speed at a stream-wise station, x , measure velocity profiles, both mean and fluctuating component, in the turbulent wake. The origin of x coordinate is at the center of the cylinder and it is increasing in the stream-wise direction.
- h) Plot velocity profiles, both mean and rms, in normalized forms.

Questions:

- 1) How does constant temperature hot-wire anemometer work?
- 2) How do you calculate mean and rms velocity components from hot wire signal?
- 3) What is the nature of the hot-wire signal in the free stream?
- 4) What do you understand by the intermittency in a turbulent wake?

Experiment # 2: Part-B

Study the dependence of shedding frequency in the Karman vortex street of the circular cylinder wake on Reynolds number.

Learning Objectives:

To understand the principles of a Constant temperature anemometer and thereby to study the vortex shedding frequency of a cylinder wake and to observe the variation of Strouhal number with Reynolds number.

Proposed Plan:

- a) You are provided with four circular cylinders of different diameters. Measure the diameters with the help of a digital Vernier.
- b) Mount the hot wire in the near wake (about 3 diameters) of the circular cylinder. Turn the wind tunnel on. At an appropriate wind speed the hot wire signal will show a sine wave. Note the frequency of the sine wave. This corresponds to the frequency of the vortex shedding.
- c) Increase the Reynolds number and note the variation in the shedding frequency.
- d) Calculate Strouhal number, $S = fd/\nu$, and Reynolds number. Plot results as $S = f(\text{Re})$.

Questions:

- 1) What do you observe when the hot wire probe is moved across the cylinder wake? Explain the reason.
- 2) Can you use this method to measure low wind speed? Explain.
- 3) Sketch the flow field on a circular cylinder as a function of Reynolds number. Start with the inviscid flow. Comment on drag and other characteristics.
- 4) What happens to the shedding frequency when the hot wire probe is moved across the cylinder near wake? Explain.

References:

1. Schlichting, "Boundary layer theory," McGraw Hill, New York, 1970.
2. Roshko, "On the development of turbulent wakes from vortex streets," NASA technical note 2913, March, 1953.
3. Turbulence - Hinze
4. A first course in turbulence - Tennekes & Lumley
5. L. Ong and J. Wallace, "The velocity field of the turbulent very near wake of a cylinder," Experiments in Fluids, 20 (1966).